


Sealed 18:1 Handpiece



**No other surgical
handpiece can provide
maintenance-free usage**

The only sealed latch head surgical implant hand piece that requires no lubrication or cleaning maintenance. Now simplicity comes along with the exceptional torque and dependability of the DynaSurg hand piece. Latch head provides internal/external irrigation.

The sealed 18:1 hand piece has the standard E-type connection and will fit most implant surgical motor systems. This hand piece is autoclavable and has a one year warranty.



*** No Maintenance Required**
*** Longer Life**

ITL DENTAL

31 Peters Canyon, Irvine, CA 92606

PH 800.277.0078 949.223.8950 FX 949.223.8960

www.itldental.com sales@itldental.com