


Bone Expander Kit

5 bone expander diameters



Ratchet Wrench



Pilot Drill



Ratchet Wrench Extender



Manual Knob



Latch Style Hand Piece Adaptor

(Kit Components Above- Not to scale)

An Alternative to Osteotomes

Bone expander drills are an alternative to osteotomes for the expansion and condensing of the atrophic mandible and maxilla in preparation for dental implant insertions. Expanders are also an alternative to the maxillary sinus elevation technique.

Surgical Trauma Decreased

Expanders are driven into the bone with a ratchet wrench or low speed hand piece. This decreases the surgical trauma of osteotomes. Bone expanders improve the clinical success by improving stability, maintaining bone density and increasing fixation.

What the Kit Includes

The Bone Expander Kit from ITL Dental includes 5, bone expanders, a ratchet wrench, a pilot drill, a ratchet wrench extender, a hand piece latch adaptor and a manual knob for finger usage. The kit is autoclavable. The cost of the kit is under \$500 making it one of the most economical expander kits.

ITL DENTAL

For more information contact: sales@itldental.com www.itldental.com or call: 800.277.0073
31 Peters Canyon, Irvine, CA 92606 949.223.8950 Fax: 949.223.8960