




Perhaps it's time to consider
the **Perio-LIG™** Syringe
from ITL DENTAL



The intraligamental syringe

PDL Injection Advantages Over Mandibular Block Injections

Patient Comfort

The cheek, tongue and lip numbness associated with conventional administration of the anesthetic is eliminated. Can be used with Vibraject™ to block the pain of dental injections.

Chair Time Efficiency

Procedure can start immediately, No 5 minute or longer waiting time.

Convenience

Precision 0.2cc dosing is provided with the syringe. Benefits are minimal anesthetic usage. Pistol style construction is durable and comfortable to use. Syringe is autoclavable.

**Perio-LIG™ Syringe for
Periodontal Intraligamental Anesthesia**

ITL DENTAL

31 Peters Canyon
Irvine, CA 92606 949.223.8950
Fax: 949.223.8960 800.277.0073
www.itldental.com sales@itldental.com