

# Gate Control Theory

Vibraject removes pain from injections based on the Gate Control Theory. The battery operated motor attaches to a conventional syringe causing the needle to vibrate. The high frequency vibration of the needle stimulates the nerve endings and blocks the transmission of pain feelings to the brain.



## Vibraject® KIT



**Includes:** ▶ One motor ▶ Three clip brackets ▶ Six batteries  
(autoclavable)

## Vibraject®

The injection syringe analgesic

- ▶ No topical anesthetic necessary
- ▶ Standard injection technique used
- ▶ Attaches to conventional syringes including intraligamental syringes
- ▶ Easy to use

Patent No: 5,647,851; 6,602,229

**ITL DENTAL**

[www.itldental.com](http://www.itldental.com)

31 Peters Canyon | Irvine California 92606 | 800-277-0073 949-223-8950 | FAX 949-223-8960 | [sales@itldental.com](mailto:sales@itldental.com)