

Oh-No



## Without **VibraShield™** to Protect



1 Assemble needle on the syringe in the normal way.

2 Slide the VibraShield™ over the needle sheath with the concave side facing the needle's point. Seat firmly against the shoulder of the sheath.

3 Remove the needle sheath with the VibraShield™ just prior to the injection.

Needle sticks are rare occurrences, but in today's world they are extremely dangerous. CDC recommends, "An engineered sharps injury protection device should be employed for recapping needles between uses and before disposal." Besides the potential of acquiring HIV/AIDS or hepatitis, the required precautionary treatment for percutaneous injuries is expensive and time consuming. For a copy of the Exposure Incident Follow-up Procedures and forms from the OSHA Bloodborne Pathogens Rule, go to: [www.osap.org/issues/pages/sharps/sharps-qna.htm#BPSareas](http://www.osap.org/issues/pages/sharps/sharps-qna.htm#BPSareas)

A Vibraject™ LLC Product

**ITL DENTAL**

31 Peters Canyon | Irvine California 92606 | 800-277-0073 949-223-8950 | FAX 949-223-8960 | [sales@itldental.com](mailto:sales@itldental.com) | [www.itldental.com](http://www.itldental.com)