


Make Your NobelBiocare™ Or ITI Straumann® Torque Wrench –Universal– for All Implant Systems.

Stop purchasing multiple torque wrenches for each implant system that you use



ITI Straumann®

ITI Straumann® Adaptor



ITL –Universal– Adaptors



NobelBiocare™ Adaptor

ITL Dental offers an adaptor that fits into the NobelBiocare™ or ITI Straumann® torque wrench that allows you to use the standard 4mm square head drivers used by most implant systems. Now, your NobelBiocare™ or ITI Straumann® torque wrench can be used with other implant systems.

If you need drivers with the 4mm square head, ITL Dental offers 43 different drivers to fit the needs of most implant systems. If you need torque wrenches, we offer the Dynatorq Torque® Wrench in settings of 10, 20, 30, 32 and 35Ncm.

Don't spend extra money when it's not necessary

To review the list of torque wrenches and drivers used by more popular implant systems, please visit the ITL Dental web site, www.itldental.com and review the chart under torque wrenches.

If you need additional assistance, please call us at: 800.277.0073

ITL DENTAL

www.itldental.com sales@itldental.com

31 Peters Canyon, Irvine, CA 92606 949.223.8950 Fx: 949.223.8960