

ITL Torque & Ratchet Wrench DRIVERS



DRIVERS

A number of drivers and accessories are available for 4mm Square Insert torque or ratchet wrenches.

DRIVERS* include HEX, TORX, SQUARE and SLOT. Drivers are available in short (posterior) and long (anterior) lengths. Drivers are available for most implant systems, making your torque or ratchet wrench universal.

EXTENSIONS The EX extension increases the driver length by 3/8 inch and the EXL by 5/8 inch.

THUMB KNOB TK This convenience knob is used with driver for manual usage. The Screw or abutment is tightened with finger pressure and final seating is accomplished with the wrench.

* For a complete sizing list of drivers, see www.itldental.com under DynaTorque® Wrench.

ITL Dental offers an array of wrench drivers for dental clinicians and implantologists

From Torx, Hex and Square drivers to Specialty drivers, our inventory is comprehensive. The source materials from which our drivers are produced and engineered, are of the finest quality.

Our various drivers are compatible with most implant systems. They are easily adapted for use on NobelBiocare™ and ITI Strumman® torque wrench products utilizing the ITL Dental adaptor. The ITL driver selection is also available in kit form and can be configured to your personal requirements.

For more specific information on our wide selection of torque wrench drivers, please contact us at any of the below-listed addresses or phone numbers.

Extensions and Accessories



EXTENSIONS



THUMB KNOBS



HEX DRIVERS



TORX DRIVERS



SQUARE DRIVERS



SLOT DRIVERS

ITL DENTAL

www.itldental.com
sales@itldental.com

31 Peters Canyon, Irvine, CA 92606 949.223.8950 800.277.0073