


# ***4mm Square Screw Driver Handle***

## ***from ITL DENTAL***



One of ITL's newest products is the 4mm Square Screw Driver Handle. This stainless steel driver handle is approximately 5.5 inches in length and accommodates the full selection of 4mm square drivers. This instrument assists in dental procedures involving the placement of implants and other related components when working in the maxilla or sinus.

Four grooves in the driver handle ensure positive-grip action when either tightening or loosening various implant components. This single unit instrument is fully autoclavable.

**ITL DENTAL**