


RIDGE SPLITTING CHISELS


ITL DENTAL is pleased to announce the following addition to our line of implantology products:

- 4 Ridge Splitting Chisels • In 2 Blade Thicknesses
- With 2mm Increment Score Lines





Ridge Splitting Chisels are available in two widths: 3mm and 6mm. Each of the two widths also are available in two blade thicknesses: 0.50mm and 2mm. Each blade is calibrated with score lines marked at 2mm increments. All chisels are available individually or as a set depending upon your specific requirements. The 4-chisel set is available with or without a fully autoclavable cassette. For product pricing contact ITL DENTAL via any of the avenues listed at the bottom of this flyer.

Score Lines



Blade Thicknesses



Ridge Splitting Chisel Autoclavable Cassette

ITL DENTAL

31 Peters Canyon, Irvine California 92606

www.itldental.com

sales@itldental.com

800.277.0073

949.223.8950

949.223.8960 Fax