


Implant Instrumentation

Implant Drivers, Extensions & Adaptors

SHAFT LENGTHS:

SHORT: 6.4mm

LONG: 12.7mm



TIPS	HEX DRIVERS	SPECS	TIPS	HEX DRIVERS	SPECS	TIPS	SPECIALTY DRIVERS	SPECS
		.035 Hex Tapered Short: HCT Long: HCLT			.070 Hex Long: HTSI			.048 / .050 Hex Short: HNHTS Long: HNHTL
		.035 Hex Short: HC Long: HCL			.080 Hex Short: HU Long: HUL			Torx Short: PROC-S Long: PROC-L NobelBiocare
		.040 Hex Tapered Short: HSM Long: HSML			2.5mm Hex Short: H25S Long: H25L			Torx Short: ITIS Long: ITI ITI/Straumann
		.048 Hex Tapered Short: HNT Long: HNLT			Hex Driver for internal Hex Long: IHDL			4mm Square Head Triangle Driver Short: TRIS
		.050 Hex Tapered Short: HTT Long: HTLT			Hex Driver for internal Hex Short: IHDS			4mm Square Head Triangle Driver Long: TRIL
		.052 Hex Tapered Short: HT2SMT Long: HT2SMLT			Implant Driver 5.7mm internal Hex Long: IHD5L			Female Driver .080 Hex Short: HFD
		.062 Hex Short: HS Long: HSL			Implant Driver 5.7mm internal Hex Short: IHD5S			TG Post for 3i/ITI/Straumann TGDN
								Mini Implant Driver Short MIDS
								Mini Implant Driver Long MIDL

(Product photos not to scale)

To order call: 800.277.0073
or contact us at: sales@itldental.com
For additional implant-related product information, visit our website at: www.itldental.com


Implant Instrumentation

Implant Drivers, Extensions & Adaptors

SHAFT LENGTHS:

SHORT: 6.4mm

LONG: 12.7mm



TIPS	LATCH STYLE DRIVERS	SPECS	TIPS	LATCH STYLE DRIVERS	SPECS	TIPS	ADAPTORS	SPECS	
		2.5mm Latch Style Hex Driver Short: H25LS Long: H25LL			Implant Latch Driver for Internal Hex Long: IHDRAL			NA NobelBiocare Torque Wrench Adaptor	
		Latch Style Torx Driver NobelBiocare Short: PROCLSS Long: PROCLSL			Implant Latch Driver for Internal Hex Short: IHDRAS			ITIA ITI/Straumann Torque Wrench Adaptor	
		.040 Latch Style, Tapered Hex Driver H40FA			Implant Latch Driver for 5.7mm Internal Hex Long: IHD5RAL			LS Latch to Square Adaptor	
		.048 Latch Style, Tapered Hex Driver T10/24			Implant Latch Driver for 5.7mm Internal Hex Short: IHD5RAS	MANUAL KNOB		SPECS	
		.050 Latch Style, Tapered Hex Driver Short: HTRAS Long: HTRAL	TIPS		SQUARE DRIVERS	SPECS		Manual Knob TK	
		4mm Latch Style Triangle Driver TRIA			.050 Square Driver, Long: CERL Short: CERS	TIPS		EXTENSIONS	SPECS
		ITI Latch Driver ITILSS Short ITILSL Long ITI/Straumann			.050 Latch Style, Square Driver T17/24			3/8" Extension EX	
			TIPS		SLOT DRIVERS	SPECS			5/8" Extension EXL
					0.017 Slot Driver Short: SN Long: SNL	<p>To order call: 800.277.0073 or contact us at: sales@itldental.com For additional implant-related product information, visit our website at: www.itldental.com</p>			

(Product photos not to scale)