


ITL DENTAL

Implant Instrumentation

Adjustable Torque Wrench



Precise Adjustable Torque Wrench with Nine Settings:
15, 20, 25, 30, 35, 40, 45, 50 and 60 Ncm
ATWP: with Precise logo ATWP: ATW: without Precise logo
4mm Square Connection

Non-Adjustable Torque Wrenches



Single Setting Torque Wrenches: 15, 20, 25, 30 and 35 Ncm
15Ncm: GTW 15; 20Ncm: GTW 20; 25Ncm: GTW 25;
30Ncm: GTW 30; 35Ncm: GTW 35
4mm Square Connection

(Above photos not to scale)

Implant Instrumentation
continued on next page...

ITL DENTAL